# Introduction to Neurofeedback Using NeuroGuide

Seamless integration of qEEG and Neurofeedback in NeuroGuide

Robert W. Thatcher, Ph.D.

# What is Neurofeedback?



- Examples:
  - Reducing Stress via Heart Rate or Blood Pressure Training
  - Achieving Calm by Galvanic Skin Response Conditioning
  - Controlling Pain by Changing Muscle Hypertension (EMG)
- <u>Neurofeedback</u> is Biofeedback of Information Contained in the Collection of Biometric Signals That Constitute the EEG, or Electroencephalogram

# Why do Neurofeedback?




- Neurofeedback is used to treat behavioral health issues by training the patient to alter their brain function via altering aspects of their EEG
- Neurofeedback has been shown to be an effective treatment modality for affective disorders (anxiety, depression, PTSD), developmental disorders (AD/HD, LDs and ASDs), acquired brain injury (CVAs and TBIs), and many other conditions related to brain function (epilepsy, chronic pain, addiction, etc.)

# What is Neurofeedback?





 Originally, Neurofeedback was done by training patients to alter aspects of their actual brainwave signals ("Raw" EEG):



# What is Neurofeedback?

 Today, Neurofeedback more commonly refers to training patients to change aspects of their quantitative EEG ("qEEG") brain maps that have been shown to be causally related to their symptoms






# What is NeuroGuide?




- An FDA Registered Database of Normative qEEG Values
- A computer program that can compare the qEEG values of a given patient to the normative database and comparative produce brain maps
- A system for collecting and recording ("acquiring") EEG in real time
- A set of tools for doing Neurofeedback on patients to help them retrain the brain functioning in order to address their behavioral health problems
- A research tool for doing statistical analyses on groups of brain maps
- Many other uses



# What is NeuroGuide?



# Why use NeuroGuide?



- Trains Only Z-scores of qEEG metrics, not absolute values
- Symptom driven doesn't "chase dysregulations"
- Relies on full 19 channel EEG
- Analyzes and rewards qEEG normalization using 3-D real time brain imaging (LORETA)
- Based on the Node-Network Hypothesis of Cortical Functioning
- Can train Surface Power, Connectivity, Current Sources (3-D) or a combination of these

# **Workshop Goals**




- To Learn and Understand a Comprehensive and Rational Approach to Analyzing and Changing Brain Function Using NeuroGuide
- To Teach Attendees How to Assess Patients, Acquire EEG an Automatic Clinical Report (ACR)
- To Learn How to Link Symptoms to Brain Dysregulations
- To Assess Symptoms By Functional Domain and Create Treatment Protocols Using NeuroLink
- To Learn How to Treat Patients Safely and Effectively as Quickly as Possible (one to two weekends) using Neuroguide and various EEG Collection Systems


# Part I: Some Theory and Background

Understanding The "Why" and the "How" of qEEG-Guided Neurofeedback Using NeuroGuide

# Linking The Patient's Symptoms to The Patient's Brain (LSB)



### Essential Steps in Helping Patients with Neurological/Psychological Problems



### **History and reason for Treatment**



### Structured Intake Interview

- What Brings You Here?
  - Symptom(s) h/x, Onset, Duration, Severity, Previous d/x, etc.
- Medical History

0

- Guided Questions re: Developmental h/x, Medical Issues, Surgeries, Physical/Psych. Trauma, TBIs, Seizures, R/Xs, Other
- Biopsychosocial Info.
- Family History and Living Situation
- Handedness, Languages
- \*\*\* Goals of Therapy \*\*\*

### Normative Testing as Deemed Relevant

#### Single Symptom Inventories / Behavior Scales

- BDI-II
- BAI
- CB-OCD
- MASC
- Etc.

Comprehensive Profile Inventories (Self and Others)

- CBRS
- PAI/MMPI
- TSI
- BPAD,
- NeuroLink\*

#### Performance-Based Screening Batteries

- CNS-VS
- NeuroTrax
- NEPSY-II
- Others

- Neuropsychological Performance Tests
  - CPT-3-CATA
  - IVA-TOVA
  - WISC/WAIS/Raven's Matrices
  - WMS-CVLT-RCFT
  - WASI/WJ-III

# Report Results and Linking Symptoms to Brain Dysregulations (LSB)

- What We're Here For:
- Assessment and NeuroLink
- Hypothesis Formation
  - EEG Collection and Recording
    - Wet vs. Dry, Cap vs. Free Sensors,
  - Recording Hygiene
    - Identifying and Eliminating Artifact: The Earlier, The Better
  - Generating The Analysis
    - Running Custom Analyses vs. Using the ACR

# Report Results and Linking Symptoms to Brain Dysregulations




Hypotheses: Linking Symptoms to Brain

Assessment Results Automatic Clin Report

### Linking Symptoms to Brain







# Analyses, Hypothesis Testing and Report Creation



# Analyses, Hypothesis Testing and Report Creation





and the local division of the local division

HERITAGE LARGE AND A

- •No Delays with Minimal Expense for a Professional Quality In-House QEEG Clinical Report
- •Less than One Minute to Produce a Professional QEEG Clinical Report, in Microsoft Word format
- ACR Provides: Empowerment, Simplicity, Accuracy & Efficiency!
- •Get Valid Normative Database Comparisons using without Depending on Internet Q-EEG Report Services!
- •Get Relevant Content and Displays, plus Helpful NFB Recommendations in Less than a Minute.
- Increased Productivity by at Least 10 Fold, e.g. Ten Reports in an Hour!

#### NeuroLink to a Clinical Report and Neurofeedback Protocol

#### **ACR Options Panel**

Only Select	Automatic Clinical Report						
Discriminant	Header	0ptions					
Functions if Patient has		— Clinic Report Options –					
a History of mTBI	Discriminant Functions						-
and/or	IBI Discriminant Function (use only if patient has a history of TBI & age > 13)						
Academic Problems	LD Discriminant Function (use only if patient has a history of academic problems & age 5 to 18)						
	Brain Performance Index (age 5 to 18)						
Select Desired		Neurofeedback Rec	ommendations				_
Neurofeedback		Addiction	Anxiety	Attention	Compulsive	Depression	
Recommendations		Secutive	Memory	Pain	Sensory	Social	
		Highest Z Score:	s				
User is Required to Certify that Only Artifact Free and High Quality Data		Z Score Threshold					
is to be Used	ertify that artifact free data was selected						
Generate Report Cancel							
				1			

Header Panel is not Active until User Certifies Data Quality


Click Generate Report after Selecting the Desired & Required Items


Network Linked to Symptoms → List of Networks appears – user selects the network, the significant surface & LORETA Z scores from the analysis are sorted and matched to the surface and Brodmann areas in the Neuroguide symptom check list and then printed in the report – if none then no NFB recommendations are produced with a statement – No Matches Present

# Part II: EEG Acquisition and Analysis


How to Acquire Good EEG, Remove Artifact, and Create qEEG Maps Using NeuroGuide


- We will:
  - Request a Volunteer to Have their EEG Recorded
  - Use NeuroGuide with a Wireless Dry Electrode System to Acquire Raw EEG
  - Use NeuroGuide to Remove Artifact and Select Brain-based EEG for Analysis
  - Use the Automated Cinical Report Generator to Create a qEEG Report





Template Artifact Rejecti	on	Z Score Artifact Rejection			
Drowsiness	High	Orowsines	s Sensitivity		
Eye Movement	Medium	Eye Mover	ment 🛓 High		
	© Low	Muscle			
Template EEG Selection		Z Score EEG S	Selection		
Amplitude Selection		Amplitude	Selection		
Amplitude Multiplier		Z Score Three			
▲ 1.00		2.00	↓ ↓		
Use Time Constraint	s Selection S	tart Time Select	tion Duration		
	≜ 00 :		3 : 🔺 00		
Generate Edits	<u>C</u> lear E	dits	Qlose		





### Automatic Clinical Report Writer Header Page

Header Options		-0	
Company Info	ABF Behavioral Health	▼	No Company Logo in Report
Company Name	ABF Behavioral Health		
Address Line 1	10840 Sheldon Road		
Address Line 1	Tuo4u Srieldon Road		
Address Line 2	Suite B		Manual Signature
			John Doe
City	Tampa State	FL	
Zin Code	33626 Country	LICA	
Zip Code	33626 Country	USA	Clinician Name (e.g. John Doe) William A. Lambos
Phone	813.235.4270 Fax	813.	
			Title (e.g. Ph.D., QEEG-D, BCIA, ECNS)
Email			Ph.D., BCN
		C 14	
	2	ave Company Info	Reset
		Exit Nex	xt


### **Automatic Clinical Report Writer Options page**

<u>H</u> eader	Options					
	— Clinic Repor	t Options				
	Discrimin	ant Functions				
	I TBI C	Discriminant Functior	n (use only if patient	t has a history of TBI &	age > 13)	
	🔘 LD D	iscriminant Function	(use only if patient	t has a history of acade	mic problems & age 5 to 1	18)
	Brain Performance Index (age 5 to 18)					
	Neurofee	dback Recommend	ations			
	Addic	tion C	Anxiety	Attention	Compulsive	Depression
	Exec	utive C	Memory	Pain	Sensory	Social
	Highest Z Scores					
	Z Score 1	hreshold				
	<b>1.96</b>					
	<b>@</b> 10	certify that artifact fro	ee data was selected	I		
			Generate	Report	ancel	

### **Automatic Clinical Report Writer Output Page**

#### Conventional EEG Samples and Quantitative EEG Analyses


Example of Linked Ears EEG and Absolute Power - Eves Closed Condition



Example of Laplacian EEG and Absolute Power - Eyes Closed Condition

00:24

FPI-CSD HUN MUMANIA and make when you and the well had a will have been and the property and the provident the pro min many many many FP2-CSD Allhach Modelat F3CSD and human man man production of the line of the production o FLCSD manufantan Anona Manufantan manuning and Manufantan manufantan C3-CSD C4-CSD P3-CSD P4-CSD O1-CSD mmemory manenon 02-CSD F7-CSD F&CSD T3-CSD T4-CSD T5-CSD proming morning T8-CSD Fz-CSD Cz-CSD warder warder warder and the second an man and a second of the second s Pz-CSD mannin 00:25 00:26 00:27 00:29 00:24 00:28



Absolute Power (uV So)

15 15 20 Frequency (Hz)


Z Score of Absolute Powe

# Part III: Linking qEEG Results to Symptoms

How to Link qEEG Mapping Results to The Patient's Symptoms

### Testing To Relate Symptoms To Patient's Brain

# NeuroLink by Applied Neuroscience, Inc.




© 2015 – All USA and International Rights Reserved


www.anineurolink.com

Press Any Key to Continue...

### NeuroLink Severity Scores By Symptom Domain




### NeuroLink and NeuroGuide Integration: Linking Symtoms to the Brain



# How Have We Learned About Brain Function?


- Open and Closed Head Wounds
  - Phinneus Gage, CVAs, TIAs, Birth Trauma, etc.
- Post-Mortem Histology Studies
  - Hippocrates, Galen, Ramon y Cajal, etc.
- Animal Ablation Research
  - Lashley, etc.
- Brain Imaging
  - Structural, Function, Axial, Proxy-based, Physiological






#### Brodmann Areas: Cytoarchitecture ("Structure") of the Brain

#### **Original Brodmann Map - Colorized**



**Outlines - with Functional Attribution** 





#### **Functional Neuroimaging Modalities**




Log\_(cm) Cerebral Cortex 2 - (tas) Lobes (occipital) 1 (10cm) **fNRs** EEG Areas (visual) EEG D (fem) **MRI** PET/SPECT Sub-areas (V1) -1 11 ...... Columns (ocular dominance) Mesoscopic Level -2 - (100,000) Minicolumns Layers -3 Cells 110...... Microscopic Neurons/Astrocytes Level Synaptosomes Molecules -5 Log-(s) ማ 5 milliseconds seconds minutes hours

#### Functional Neuroimaging Modalities


# Traditional Domains of Neuropsychological Functioning and Assessment



# Six Functional Modules as Measured by fMRI




#### Six Functional Modules as Measured by PET/SPECT




#### Six Functional Modules as Measured by qEEG LORETA


#### Hagmann et al. Modules




MOD 3









MOD 5






#### LORETA Effective Connectivity: Magnitude and Direction of Information Flow Between Brodmann Areas







# Information Flow in the Scalp Surface EEG



#### Information Flow: The Node-Network Hypothesis of Brain Function




# The Process of Neuropsychological Assessment





# **Electrical Neuroimaging & Cortical Source Localization**

#### **Cortical Surface Projection**

#### Horizontal, Sagital & Coronal Views of a Single Slice









# Electrical Neuroimaging – Assessment and Treatment


- Treating the patient:
  - Preparing for 19-channel Z-scored Neurofeedback
  - Loading the Protocol
  - Choosing Surface vs. LORETA NFB
  - Train Symptoms or Networks?
  - NFB vs. BrainSurfer
  - Choosing Parameters:
    - Feedback Displays
    - Thresholds
    - Durations
  - Monitoring Progress
    - Within the Session
    - Between Sessions
    - Re-checking symptomatology with NeuroLink history









		1 1	
	Network	Function	Disorder
	Hippocampal- diencephalic and parahippocampal- retrosplenial	<ul> <li>memory</li> <li>spatial orientation</li> </ul>	•Amnesias •Korsakoff's syndrome •Mild Cognitive impairment •Alzheimer's disease (early) •Balint syndrome
	Temporo- amydgala- orbitofrontal	<ul> <li>Behavioural inhibition</li> <li>Memory for temporally complex visual information</li> <li>Olfactory-gustatory-visceral functions</li> <li>Multimodal sensory integration</li> <li>Object-reward association learning</li> <li>Outcome monitoring</li> </ul>	<ul> <li>Alzheimer's Disease (advanced)</li> <li>Semantic dementia</li> <li>Klüver-Bucy syndrome</li> <li>Temporal lobe epilepsy</li> <li>Geschwind's syndromes</li> <li>Psychopathy</li> <li>Bipolar affective disorders</li> </ul>
encephalic and al-retrosplenial	Dorsomedial default network	Pain perception     Self-knowledge     Attention     Mentalizing	Depression     Autism     Schizophrenia     Obsessive compulsive disorder     Mild Cognisius Imposisments
a- ork work'		•Empathy •Response selection and action monitoring •Autobiographical memory •Person perception	<ul> <li>Mild Cognitive Impairmentt</li> <li>Alzheimer's Disease (early)</li> <li>Attention Deficit Hyperactivity Disorder</li> <li>Anxiety</li> </ul>







### **LORETA** Coherence



## **LORETA Absolute Phase**



### **LORETA Magnitude and Direction of Information Flow**




# Loreta Default Brain Network




#### 


#### Temporal




#### Posterior Cingulate




#### Anterior Cingulate



#### Hippocampus



#### Parietal



# **Some Relevant Items and Questions**

- Loops in the Brain and Why Homeostasis and Equilibrium are Critical for Brain Function
- The EEG is Produced Exclusively by Summated Synaptic Potentials
- How does EEG Biofeedback Change the EEG?
- How does EEG Biofeedback Change Synaptic Potentials?
- What are the Mechanisms of Modification of Synapses by Operant Conditioning of EEG at the Molecular Level? (nu. Accumbens & Reinforcement


# **Essentials of Operant Conditioning**

- 1. There must be a 'real' and 'valid' neural event to be reinforced
- 2. The 'Reinforcement' must be distinct and clear
- 3. The interval of time between the spontaneous 'emitted event' and the 'reinforcement' can not be too short, approximately < 250 msec? or too long approximately > 20 sec
- 4. The Schedule of Reinforcement is Important with two General Types 'Continuous' vs 'Partial' Reinforcement -Continuous is good at the beginning but not as resistant to extinction as is Partial Reinforcement

#### **EEG Neurofeedback** Memory trace Memory trace Modulated decay of the memory trace Simple decay of the memory trace Reward Reward Firing rate Firing rate No reward 1 No reward Time Time Trial 1 Trial 2 Trial 3 Trial 1 Trial 2 Trial 3 D1 or D2 Cerebral Cortex (FEF, SEF, DLPF, LIP) antagonist Basal Ganglia CD D2 D1 GPe CD SC SNc SNr SNc STN SNr 10 SC Saccade


### Z Score NFB – Simplicity, Accuracy & Standardization

#### Difference Between Standard Neurofeedback vs Z Score Neurofeedback




Move Z toward 0






18-4 Changes underlying short- and long-term memory in a single sensory and motor neuron.











\* Combine with Neurofield in cases to reset or "unstick" the brain

### Seamless QEEG and Neurofeedback Approximately 50 – 60 minutes for a Single Session in Three Steps from Clinical Interview to QEEG to Neurotherapy



### Seamless QEEG and Neurofeedback Approximately 50 – 60 minutes for a Single Session in Three Steps from Clinical Interview to QEEG to Neurotherapy



# Select a Network or Symptoms, Frequency and Metric


Symptoms ICN	Networks	Neuropsychological	D	oD/VA		
	Network			Seve	rity	
Addiction				0		
Anxiety				0		
Attention - Dorsal				0		
Attention - Ventral				0		
Attention - Emotional				0		
Default Mode				0		
Executive Function				0		
Face, Object Recogni	tion			0		
Language				0		
Memory - Emotion				0		
Mirror Neuron				0	)	
Mood				0	)	
Pain				0	)	
Pleasure				0	)	
Salience				0	)	
Schizophrenia				0	)	
Working Memory				0	)	
DTI - Frontal Limbic				0		
DTI - Frontal Occipital				0		
DTI - Frontal Parietal						
DTI - Frontal Tempora				0		
DTI - Local Frontal				0		
DTI - Local Limbic				0		
DTI - Local Occipital				0		
DTI - Local Parietal				0		
DTI - Local Temporal				0		
Hagmann Module 1 (\	/ision)			0	)	
Hagmann Module 2 (A	ttention, Wor	king Memory)		0	)	
Hagmann Module 3 (A	uditory, Lang	uage, Memory)		0	)	
Hagmann Module 4 (A	Auditory, Lang	uage, Memory)		0	)	
Hagmann Module 5 (Executive, Sequential Planning)					)	
Hagmann Module 6 (Executive, Social Skills)					)	
Isocortex Hippocampo	Isocortex Hippocampocentric					
Isocortex Olfactocentric				0		
Mesocortex Hippocampocentric				0		
Mesocortex Olfactocentric				0		
Mesulam - Emotional I	Mesulam - Emotional Memory					
Mesulam - Executive	0					
Mesulam - Face/Object Indentification						
Mesulam - Language						
Mesulam - Spatial Atte	Mesulam - Spatial Attention					

Symptoms	ICN	Networks	Neuropsychological	D	oD/VA		
Intrinsic Connectivity Network						rity	٨
ICN 1 (Limbic	, Medial-	Temporal, Em	otion)		0		
ICN 2 (Rewar	rd, Emoti	on)			0		
ICN 4 (Langu	age, Exe	cutive)			0		
ICN 6 (Premo	tor, Supp	lemental Mot	or)		0		
ICN 7 (Visual-Spatial Processing)					0		
ICN 8, 17 (Primary Sensory Motor)				0			
ICN 9 (Parietal)					0		
ICN 10 (Picture Naming, Visual Tracking)					0		
ICN 11, 12 (V	lisual Sys	tem)			0		
ICN 13 (Defa	ult Mode	Network)			0		¥
ICN 15 (Right Hemisphere, Attention, Reasoning, Memory)				0			
ICN 16 (Auditory, Music)					0		
ICN 18 (Left Hemisphere, Language)					0		v
	ICN 1 (Limbic ICN 2 (Rewar ICN 4 (Langu ICN 6 (Premo ICN 7 (Visual ICN 8, 17 (Pri ICN 9 (Pariett ICN 10 (Pictu ICN 10 (Pictu ICN 11, 12 (V ICN 13 (Defa ICN 15 (Right ICN 16 (Audit	Intrins ICN 1 (Limbic, Medial- ICN 2 (Reward, Emotia ICN 4 (Language, Exe ICN 6 (Premotor, Supp ICN 7 (Visual-Spatial F ICN 8, 17 (Primary Ser ICN 9 (Parietal) ICN 10 (Picture Namin ICN 11, 12 (Visual Sys ICN 13 (Default Mode ICN 15 (Right Hemisp) ICN 16 (Auditory, Mus	Intrinsic Connectivi ICN 1 (Limbic, Medial-Temporal, Em ICN 2 (Reward, Emotion) ICN 4 (Language, Executive) ICN 6 (Premotor, Supplemental Mot ICN 7 (Visual-Spatial Processing) ICN 8, 17 (Primary Sensory Motor) ICN 9 (Parietal) ICN 10 (Picture Naming, Visual Trac ICN 10 (Picture Naming, Visual Trac ICN 11, 12 (Visual System) ICN 13 (Default Mode Network) ICN 15 (Right Hemisphere, Attentio ICN 16 (Auditory, Music)	Intrinsic Connectivity Network ICN 1 (Limbic, Medial-Temporal, Emotion) ICN 2 (Reward, Emotion) ICN 4 (Language, Executive) ICN 6 (Premotor, Supplemental Motor) ICN 7 (Visual-Spatial Processing) ICN 8, 17 (Primary Sensory Motor) ICN 9 (Parietal) ICN 10 (Picture Naming, Visual Tracking) ICN 11, 12 (Visual System) ICN 13 (Default Mode Network) ICN 15 (Right Hemisphere, Attention, Reasoning, Memory) ICN 16 (Auditory, Music)	Intrinsic Connectivity Network ICN 1 (Limbic, Medial-Temporal, Emotion) ICN 2 (Reward, Emotion) ICN 4 (Language, Executive) ICN 6 (Premotor, Supplemental Motor) ICN 7 (Visual-Spatial Processing) ICN 8, 17 (Primary Sensory Motor) ICN 9 (Parietal) ICN 10 (Picture Naming, Visual Tracking) ICN 10 (Picture Naming, Visual Tracking) ICN 11 (Default Mode Network) ICN 15 (Right Hemisphere, Attention, Reasoning, Memory) ICN 16 (Auditory, Music)	Intrinsic Connectivity Network         Seve           ICN 1 (Limbic, Medial-Temporal, Emotion)         0           ICN 2 (Reward, Emotion)         0           ICN 4 (Language, Executive)         0           ICN 6 (Premotor, Supplemental Motor)         0           ICN 7 (Visual-Spatial Processing)         0           ICN 8, 17 (Primary Sensory Motor)         0           ICN 9 (Parietal)         0           ICN 10 (Picture Naming, Visual Tracking)         0           ICN 13 (Default Mode Network)         0           ICN 15 (Right Hemisphere, Attention, Reasoning, Memory)         0           ICN 16 (Auditory, Music)         0	Intrinsic Connectivity Network       Severity         Intrinsic Connectivity Network       Severity         ICN 1 (Limbic, Medial-Temporal, Emotion)       0         ICN 2 (Reward, Emotion)       0         ICN 4 (Language, Executive)       0         ICN 6 (Premotor, Supplemental Motor)       0         ICN 7 (Visual-Spatial Processing)       0         ICN 8, 17 (Primary Sensory Motor)       0         ICN 9 (Parietal)       0         ICN 10 (Picture Naming, Visual Tracking)       0         ICN 13 (Default Mode Network)       0         ICN 15 (Right Hemisphere, Attention, Reasoning, Memory)       0         ICN 16 (Auditory, Music)       0

Symptoms	ICN	Networks	Neuropsychological	D	oD/VA	
Neuropsychological Diagnosis						ity <sup>′</sup>
Agnosia of Ac	ction App	perceptive			0	
Agnosia of Ac	ction Ass	ociative			0	
Agnosia Audi	tory App	erceptive			0	
Agnosia Auditory Associative					0	
Agnosia Auditory Space					0	
Agnosia Prosopagnosia (Face)					0	
Agnosia Social Emotional				0		
Agnosia Social of Action - Theory of Mind				0		
Agnosia Somatosensory Autotopagnosia				0		
Agnosia Somatosensory Finger				0		

Symptoms	ICN	Networks	Neuropsychological	D	oD/VA	
	Seve	rity				
Attention - Re	-Experie	nces Intrusive	e Memories		0	
Attention - Em	notional	Numbing			0	
Attention - Dis	stracting	Pain			0	
Attention - Difficulty Multi-Tasking						
Attention - Worsens with Emotional Stress						
Attention - Dissociative Episodes				0		
Attention - Worsens With Withdrawl Symptoms				0		
Chronic Pain - Neuropathic				0		
Chronic Pain - Musculoskeletal				0		
Chronic Pain - Diffuse Pain (Entire Body)				0		


Symptoms ICN Networks Neuropsychological I	DoD/VA
Symptom / Complaint	Severity
nosognosia - Denial of a Problem	0
nxiety	0
ttention Deficits - Easily Distractible	0
uditory Sequencing Problems	0
alance Problems	0
lurred Vision	0
hronic Pain	0
ompulsive Behaviors and/or Thoughts	0
oncentration Problems	0
ecreased Tactile or Skin Sensitivity	0
elusional	0
epression (Sad & Blue)	0
ifficulty Comprehending Social Cues	0
yscalcula - Problems Calculating	0
yslexia - Letter Reversal	0
kecutive Function Problems	0
ace Recognition Problems	0
ailure to Initiate Actions	0
	0
yperactive and/or Agitation	0
ipulsive Behaviors	0
sensitive to Others Emotional Expressions	0
sensitive to Other's Feelings	-
w Motivation	0
ow Threshold for Anger & Loss of Control	0
igrane Headaches	0
lood Swings	0
ulti-Tasking Problems	0
bsessive Thoughts about Self	0
bsessive Thoughts and/or Hyper Focused	0
ppositional Defiant Conduct	0
rientation in Space Problems	0
erception of Letters Problems	0
oor Judgement	0
oor Skilled Motor Movements	0
oor Social Skills	0
eceptive Language Problems	0
ecognizing Objects by Touch Problems	0
elf-Esteem Problems	0
equential Planning Problems	0
hort-Term Memory Problems	0
low Reader	0
lowness of Thought - Easily Confused	0
patial Perception Problems	0
peech Articulation Problems	0
peccri Auculau011 F10Diettis	0
ubstance Abuse	
ubstance Abuse ymptoms of Fibromyalgia	0



# **Neuroimaging Neurofeedback Symptom Check List**

**Click Symptoms or Neuropsychological Diagnoses** 

List of Matching




# **Use the Progress Chart as a Feedback Display**


#### Neurofeedback Setup Panel

Move the Display to the Client's Monitor


<i>ŵ</i>	Surface Neurofeedback – 🗆 🗙	
Protocol Session Rounds Progress		neurofeedback Display - 🗆 🗙
Metric       Frequency <ul> <li>Absolute Power</li> <li>Relative Power</li> <li>Power Ratio</li> <li>Amplitude Asymmetry</li> <li>Coherence</li> <li>Absolute Phase</li> <li>Phase Shift</li> <li>Phase Lock</li> </ul> <ul> <li>Montage Reference</li> <li>Linked Ears</li> <li>Average Reference</li> <li>Laplacian</li> </ul> <ul> <li>Frequency</li> <li>Delta</li> <li>Theta</li> <li>Alpha</li> <li>Beta</li> <li>Alpha 1</li> <li>Alpha 2</li> <li>Beta 1</li> <li>Beta 2</li> <li>Beta 3</li> </ul>	Surface Neurofeedback         –         ×           D/T         D/T         Auto Spectral Channels - Absolute Power	Neurofeedback Display
2.00 88 0.25 sec Progr	ess Charts Sound Fz Cz Pz Load	





#### Surface Neurofeedback Inter-Session Progress



X









# Examples of Surface EEG Changes After EEG Neurofeedback


#### **Pre-Treatment**


### Post – 10 Treatments











# **Examples of Electrical Neuroimaging After Neurofeedback**

